

Toolkit for Establishing Laws to Control the Use of Lead in Paint

Module C.i.

Analytical Methods for Measuring Lead in Blood



Outline

- Background
- Essentials of sample collection
- Brief information on different analytical methods
- Quality control considerations
- Summary
- References
- Disclaimer
- Point of Contact



Background

- Lead exposure is primarily assessed through measurement in whole blood.
- The most common laboratory methods to measure blood lead concentrations are:
 - Anodic Stripping Voltammetry (ASV)
 - Atomic Absorption Spectrometry (AAS)
 - Inductively Coupled Plasma Mass Spectrometry (ICP-MS)

Analytical methods differ in their limit of detection, accuracy, costs and technical requirements (e.g. sample preparation, calibration, and skilled personnel)



Sample collection

Care is needed

- Essential to avoid external contamination of the sample.
 - Personnel should be trained in good sampling and handling techniques to avoid contamination.
 - Collect, store and transport samples in a lead-free environment.
 - Thoroughly cleanse the skin around the puncture site.
 - Use lead-free sampling equipment and tubes. If not available send 'blanks' from same batch to the laboratory for testing of background lead content.
- Observe universal biosafety precautions.



Sample collection

Care is needed

- Collect whole blood in a tube containing EDTA or heparin.
 - Invert the filled tube 8–10 times to ensure adequate mixing.
 - Clotted samples should be rejected – analytical results will be unreliable.
- Make sure to label the tube with the patient's identification details.
- Refrigerate samples ($<4^{\circ}\text{C}$) that are awaiting analysis – do not freeze.
 - Note: does not apply to samples measured using point-of-care device, which should be kept at room temperature.



Choice of analytical method is determined by resources and needs

- Resource issues include:
 - availability of trained laboratory staff;
 - cost of reagents and other materials e.g. special gases, compressed air;
 - typical number of analyses needed (cost per analysis)
 - economy of scale possible with methods that allow multiple analyses;
 - special operating requirements e.g. reliable electricity supply, cooling water.



Choice of analytical method is determined by resources and needs

- Required limit of detection varies according to the reason for the analysis.
- Population studies – may need a method accurate to $<1 \mu\text{g/dL}$
 - e.g. geometric mean blood lead concentration in USA in 2011–12 was $0.973 \mu\text{g/dL}$.
- Confirmation of lead exposure and decisions on management – method accurate to $5 \mu\text{g/dL}$ acceptable.
 - Note: method may need to go to $>65 \mu\text{g/dL}$ in severe cases of poisoning



Anodic stripping voltammetry (ASV)

- Both laboratory-based and point-of-care devices available.
- EDTA is the preferred anticoagulant for laboratory method.
- Can analyse small samples: 50–100 μL .



Anodic stripping voltammetry (ASV)

Laboratory method

- Relatively low-cost.
- Requires skilled laboratory technician and good quality reagents for best results.
- Sample pre-treatment is needed.
- Typical analytical range is 1–100 $\mu\text{g/dL}$, but greatest precision at blood lead concentrations $>10 \mu\text{g/dL}$.
- May be interference from elevated blood copper.
- Largely superseded by other methods.



ASV: Point-of-care device

Considerations & limitations

- Portable device, can run on batteries – can be taken to the site.
- Uses a finger-prick (capillary sample), though venous samples can also be used.
- Equipment is supplied with calibration device and controls for high and low blood lead concentrations.




ASV: Point-of-care device

Considerations & limitations


- Laboratory technician is not required to perform measurement – any scientifically competent person can be trained to use the equipment.
- Result available within minutes so immediate decisions can be made about management.

(Pictures adapted from reference 1)


Collect capillary or venous sample



Put blood into a treatment reagent tube and mix by inverting 8–10 times



Place a drop of sample on sensor. Results in 3 minutes



ASV: Point-of-care device

Considerations & limitations

- Widely-used device, LeadCare II, has analytical range of 3.3 - 65 $\mu\text{g}/\text{dL}$.
- Has comparable accuracy with laboratory-based methods.
- Elevated blood lead concentrations should, however, be confirmed with a laboratory-based method.
- Some experience of using LeadCare II to measure higher blood lead concentrations by diluting the sample.
(Reference 2)



ASV: Point-of-care device


Considerations & limitations

- Risk of sample contamination is high – extra care is needed:
 - finger-prick site likely to be highly contaminated with lead and needs thorough cleansing;
 - location of exposure likely to be highly contaminated e.g. with dust, so samples should be taken and analysed in a clean room.



Atomic Absorption Spectrometry (AAS)

- Two techniques: Flame Atomic Absorption Spectrometry (FAAS) and Graphite Furnace Atomic Absorption Spectrometry (GFAAS).
- Methods differ in sample size needed, limits of detection and complexity of sample preparation.



Flame Atomic Absorption Spectrometer



Graphite Furnace Atomic Absorption Spectrometer

Flame Atomic Absorption Spectrometry (FAAS)

- Relatively easy to use and moderate cost.
- Needs special gases.
- Can be fitted with autosampler so multiple samples can be processed.
- Limit of detection depends on sample preparation and method used
 - at best: $\sim 10 \mu\text{g/dL}$ with sample size of $50\text{--}100 \mu\text{L}$.



Graphite Furnace Atomic Absorption Spectrometry (GFAAS)

- Requires skilled laboratory technician.
- Needs special gases.
- Can analyse very small samples: 10-50 μL .
- Methods available that can measure lead concentrations $<0.1 \mu\text{g/dL}$, though in routine use limit of detection is around 1–2 $\mu\text{g/dL}$.
- Can be fitted with autosampler so large number of samples can be run.
- Can be set up to measure multiple trace elements.



LEAD PAINT ALLIANCE



World Health
Organization

Inductively-coupled plasma mass spectrometry (ICP-MS)

- Expensive and has high running costs
 - more economical if used for large sample runs.
- Requires highly-skilled laboratory technician.
- Very low limit of detection: 0.1 $\mu\text{g}/\text{dL}$.
- Can measure multiple elements from a small sample (50–100 μL).
- Can determine isotope ratio, which may help to identify the source of the lead.



LEAD PAINT ALLIANCE



World Health
Organization

Lead isotope ratios can help to identify source

- Four main isotopes of lead are 208, 206, 207, 204.
- Ratios of the isotopes vary by source of the ore.
- Isotope ratio of soils represents mixing of lead from the various ores used in gasoline, consumer products and smelting.
- If isotope ratio in a lead source and in blood can be characterized, then this can be useful 'fingerprinting' of environmental pollution.

(Reference 3)



**World Health
Organization**




Chart shows group of children exposed to same source of lead and an individual exposed to a different source.

(Reference 4)



LEAD PAINT ALLIANCE

Quality control considerations

- Important that analytical results are reliable.
- Laboratory should have in place an adequate quality management system e.g.:
 - standard operating procedures;
 - documented training and monitoring of staff performance;
 - use of certified reference standards;
 - internal quality control procedures – daily checks of analytical accuracy;
 - participation in external quality assessment scheme e.g. US LAMP.



Laboratory quality assurance - LAMP

- A voluntary program that focuses on assuring the quality of blood lead, cadmium, and mercury analyses.
- Each quarter US CDC provides spiked blood samples, which are analyzed by participating laboratories who return the results to CDC.
- CDC provides detailed reports on the laboratories about how well they performed these analyses.
- No charge for participation.

The logo for the Laboratory Quality Assurance Program (LAMP) features the word "LAMP" in large, white, sans-serif capital letters. The text is set against a background of overlapping, semi-transparent yellow and orange circular shapes that create a sense of depth and movement.

Centers for Disease Control and Prevention (CDC)
Lead and Multi-Element Proficiency
4770 Buford Highway N.E., Mailstop F-18
Atlanta, GA 30341-3724 USA

Fax number: (770) 488-4097
E-mail address: LAMP@cdc.gov

Summary

- Whole blood is the preferred sample for assessing exposure to lead.
- Adequate measures should be taken to avoid sample contamination.
- A range of analytical methods are available – the decision about which one to use is determined by the available resources and the limit of detection required.
- Quality assurance procedures are important to ensure the reliability of analytical results.



References

1. LeadCare® II Blood Lead Analyzer User's Guide. N. Billerica (MA): Magellan Diagnostics, Inc.; 2015
([http://www.leadcare2.com/getmedia/73ac501b-35e3-4d74-b6a9-9e7fc51adef0/70-6551_Rev_07_User-s_Guide,_LeadCare_II_\(PRINT\)-1.pdf.aspx](http://www.leadcare2.com/getmedia/73ac501b-35e3-4d74-b6a9-9e7fc51adef0/70-6551_Rev_07_User-s_Guide,_LeadCare_II_(PRINT)-1.pdf.aspx), accessed 20 Jan 2017).
2. Neri AJ et al. Analysis of a novel field dilution method for testing samples that exceed the analytic range of point-of-care blood lead analyzers. *Int J Environ Health Res.* 2014; 24(5):418-428)
3. Komárek M et al. Lead isotopes in environmental sciences: A review. *Environment International.* 2008; 34:562–577
4. Brown MJB, US Centers for Disease Control and Prevention, personal communication, 2015



Additional references

Sample collection

Step-by-step guide for collecting capillary sample.
US Centers for Disease Control and Prevention.

Poster:

http://www.cdc.gov/labstandards/pdf/vitaleqa/Poster_CapillaryBlood.pdf

Video demonstration:

http://www.cdc.gov/nceh/lead/training/blood_lead_samples.htm

Guidelines on drawing blood: best practices in phlebotomy. Geneva: World Health Organization; 2010

(http://www.who.int/injection_safety/sign/drawing_blood_best/en/, accessed 3 January 2017).




Additional references

Analysis

Brief guide to analytical methods for measuring lead in blood. Geneva: World Health Organization; 2011 (available in Chinese, English, French and Spanish) (http://www.who.int/ipcs/assessment/public_health/lead/en/, accessed 9 February 2017)

Lead and Multi-element Proficiency programme (LAMP), [website]. Atlanta (GA): US Centers for Disease Control and Prevention (<http://www.cdc.gov/labstandards/lamp.html>, accessed 9 February 2017)



Disclaimer

© WHO 2017

Some rights reserved. This work is available under the CC BY-NC-SA 3.0 IGO licence

The mention of specific companies or of certain manufacturers' products does not imply that they are endorsed or recommended by the World Health Organization in preference to others of a similar nature that are not mentioned. Errors and omissions excepted, the names of proprietary products are distinguished by initial capital letters.

All reasonable precautions have been taken by the World Health Organization to verify the information contained in this publication. However, the published material is being distributed without warranty of any kind, either expressed or implied. The responsibility for the interpretation and use of the material lies with the reader. In no event shall the World Health Organization be liable for damages arising from its use.



**World Health
Organization**



LEAD PAINT ALLIANCE

Point of Contact

Joanna Tempowski

Department of Public Health, Environmental and Social

Determinants of Health

World Health Organization

20 Avenue Appia

211-Geneva 27

Switzerland

Email: tempowskij@who.int

Date: February 2017



World Health
Organization



LEAD PAINT ALLIANCE