





Introduction to the training movie: Passive sampling of ambient air - pre-cleaning of PUF disks

RESEARCH CENTRE FOR TOXIC COMPOUNDS IN THE ENVIRONMENT (RECETOX)

Masaryk University, Brno, Czech Republic

Stockholm Convention Regional Centre in the Czech Republic (SCRC)

Kateřina Šebková, Director of the SCRC, sebkova@recetox.muni.cz



Stockholm Convention Regional Centre for Capacity Building and the Transfer of Technology

Final Results Workshop of the UNEP/GEF project on 'Establishing the Tools and Methods to Include the Nine New POPs into Global Monitoring Plan' and Inception Workshop for UNEP/GEF project on 'Continuing Regional Support for the POPs GMP under Stockholm Convention in the Africa Region", Ghana, Accra, 6-8 July 2016

Outline



- 1. About us Who is RECETOX and what do we do?
- 2. Introduction to the training movie about passive sampling of ambient air
- 3. Movie screening "instructive movie for the cleaning of PUF disks"
- 4. Questions and Answers







RECETOX in brief



Research Centre for Toxic Compounds in the Environment

RECETOX , Faculty of Science, Masaryk University, Brno, Czech Republic


www.recetox.muni.cz

Functions

- Research and educational institution
- International Research Infrastructure trans-national open-access
- Cooperation and Capacity Building
 - √ National Centre for Toxic Compounds
 - ✓ Stockholm Convention Regional Centre for Capacity Building and the Transfer of Technology in the Czech Republic supporting in particular Central and Eastern European Region

Staff and equipment

+200 employees, 4000 m² of research facilities, 85 Ph.D. students, 70 laboratories, 13 working groups in 4 Research Divisions, state of the art facilities, in the material and other equipment, more than 30 years of expertise



About RECETOX - Structure

4 Research Divisions

- 1. Environmental Chemistry and Modelling
- 2. Organic Photochemistry and Supramolecular Chemistry
- 3. Protein Engineering
- 4. Ecotoxicology

3 Core Facilities - RESEARCH INFRASTRUCTURE

- 1. Trace Analytical Laboratories
- 2. GENASIS information system
- 3. ELSPAC database



Two internal bodies facilitating collaboration with stakeholders and marketing capacity, research outcomes and use of infrastructure:

- 1. National Centre for Toxic Compounds = technical support in implementation of chemical conventions (Basel, Rotterdam and Stockholm), SAICM and Minamata Convention on Mercury
- 2. Stockholm Convention Regional Centre for Capacity Building and the Transfer of Technology (SCRC Czech Republic)







Our work

Environmental policy, international conventions

Linking environment and health

Data evaluation, analysis, interpretation and modeling

Environmental management

Environmental fate studies

Objectives: delivery of complex information to protect humans and ecosystems from environmental hazards and linking environment and health

Information management and development of databases

Ecological and human risk assessment

Monitoring of environmental levels

Production and application of toxic chemicals

Assessment of harmful effects







Stockholm Convention Regional Centre


- Hosted by RECETOX
- established in 2007,
- endorsed in 2009, prolonged 2013, 2015
- provides support to the Stockholm Convention on POPs contracting Parties (countries) in its geographical region of the Central and Eastern Europe and beyond.
- chemicals management and management of wastes containing toxic chemicals by providing training, capacity building, expertise support in a number of fields

www.recetox.muni.cz/RC

Work also covers additional international instruments related to chemicals management – i.e. SAICM, Basel, Rotterdam and Stockholm Conventions), and Minamata Convention on Mercury



Stockholm Convention Regional Centre for Capacity Building and the Transfer of Technology











Monitoring networks



www.monet.recetox.muni.cz

MONET, MOnitoring NETwork, is a monitoring program operated by the Research Centre for Toxic Compounds in the Environment (RECETOX) of the Masaryk University in Brno, Czech Republic.

The programme aims at detection of environmental contaminants (toxic chemicals) in air, water, soil, and biota by passive sampling and spreads over the three continents - Europe, Africa and Asia.

Ambient air sampling is done by samplers with PUF disks.












Stockholm Convention
Regional Centre for Capacity Building
and the Transfer of Technology

MONET sampling sites worldwide





Chemicals covered in MONET

The MONET programme targets persistent organic pollutants, polyaromatic hydrocarbons, heavy metals, some endocrine disrupting chemicals and other toxic compounds including currently used pesticides), however the analyses differ between the individual MONET networks.

In relation to the Stockholm Convention, the MONET Africa follows a full set of the Stockholm Convention chemicals + some candidate chemicals.

All generated information is stored onlin in the GENASIS, environment data repository and portal at www.genasis.cz









Contacts for MONET activities

For more information on RECETOX's MONET activities or cooperation please contact:

Dr. Kateřina Šebková, PhD. (<u>sebkova@recetox.muni.cz</u>) - director of the Stockholm Convention Regional Centre for establishing a new cooperation

Dr. Petra Přibylová, PhD. (<u>pribylova@recetox.muni.cz</u>) - head of the Trace Analytical Laboratories for coordination of ongoing sampling and analyses of samples,

Dr. Roman Prokeš, PhD (<u>prokes@recetox.muni.cz</u>) - on sampling procedures, dispatching of materials and delivery of samples.









Introduction to the training video









Introduction to the training video

General information

The video was produced by SCRC RECETOX for UNEP Chemicals and Waste Branch in spring 2015 as a training movie with instructions for the cleaning of polyurethane foam (PUF) disks used in passive samplers for ambient air.


The work was undertaken as a part of the UNEP/GEF project 4B97 "Establishing the Tools and Methods to Include Nine New POPs into the Global Monitoring Plan" for the component on POPs analysis and capacity building in developing countries.


10 minutes long, voice-over in four languages (EN, FR, RU and ES).

Acknowledgements of involvement of the RECETOX research infrastructure laboratories and staff as well as comments received by GMP partners from Canada, Netherlands, Spain and Sweden.



Sampling device containing PUF disk - AIR











PUF = polyurethane foam

Introduction to the training video

The chapters:

- 1. Context and available guidance for monitoring under the Stockholm Convention
- Availability of the polyurethane foam (PUF) and its shaping into PUF disks
- 3. Pre-cleaning procedures of the PUF disks
- 4. Storage of the clean PUF disks
- 5. Work in the field at the sampling site (preparation, necessary material, sampling protocol)
- 6. First deployment of the PUF disk (exchange of the disk) at the sampling site and recording of supplementary information
- 7. Transport of PUF disk to laboratory for analyses









Training video - PUF disks cleaning in detail

PUF disks pre-cleaning consists of these **steps** in the laboratory:

- 1. Cutting PUF disks from a purchased material (PUF sheets)
- 2. Setting up an extraction apparatus (Soxhlet extractor) by keeping sound laboratory practice and safety rules at work
- 3. Performing the extraction and exchange of solvents
- 4. Drying the pre-cleaned PUF disks
- 5. Packaging PUF disks for controlled storage (or transport)









Screening of the training video









Are there any Questions?









Thank you for the opportunity to present our work and for your kind attention











www.recetox.muni.cz/RC

To know more about RECETOX, please visit:

www.recetox.muni.cz

VIDEO: http://youtu.be/JxmryXtdnKo

NEWSLETTER (cz-en-ru): http://www.recetox.muni.cz/index-en.php?pg=information-materials--recetox-newsletter

INFORMATION MATERIALS: http://www.recetox.muni.cz/index-en.php?pg=information-materials

www.genasis.cz www.pops-gmp.org www.elspac.cz













